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
Tangential movement with a high-speed saline jet

**Selects**

Target necrotic tissue and debris using the localized vacuum<sup>1,3-5</sup>

**Excises**

Ablate nonviable tissue with maximum precision<sup>1,3-5</sup>

**Evacuates**

Remove debris and slough while preserving viable tissue<sup>1,3-5</sup>

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


\* Saving in the average number of procedures with VERSAJET (1.91-1.26=0.65) \* average cost per debridement procedure reported in Granick, et al [ref 3] of \$3,393 in 2002/03, of which the average cost of disposables per procedure was \$802 (24%). Savings are shown before the cost of the VERSAJET hand piece

\*\* Based on five clinical studies involving 451 patients with mixed wounds (acute, chronic, trauma, burns) debrided with VERSAJET: (Granick,[ref 3]; Gravante,[ref 9]; Gurunloughlu, [ref 10]; Mosti & Mattaliano,[ref 4]; Vanwijck,[ref 6])

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
\*\* Vanwijck [ref 5]

# Références

- 1. Granick MS, Posnett J, Jacoby M, Noruthun S, Ganchi PA, Datiashvili RO. Efficacy and cost-effectiveness of a high-powered parallel waterjet for wound debridement. *Wound Repair Regen.* 2006;14:394-397.
- 2. Mosti G, Maltaliano V. The debridement of chronic leg ulcers by means of a new, fluidjet-based device. *Wounds.* 2006;18:227-237.
- 3. Granick M, Boykin J, Gamelli R, Schultz G, Tenenhaus M. Toward a common language: surgical wound bed preparation and debridement. *Wound Repair Regen.* 2006;14:S1-S10.
- 4. Cubison TC, Pape SA, Jeffery SL. Dermal preservation using the Versajet hydrosurgery system for debridement of paediatric burns. *Burns.* 2006;32:714-720.
- 5. Vanwijck R, Kaba L, Boland S, Gonzales y Azero M, Delange A, Tourbach S. Immediate skin grafting of sub-acute and chronic wounds debrided by hydrosurgery. *J Plast Reconstr Aesthet Surg.* 2010;63:544-549.
- 6. Abernathie B, Granick MS. Centrifugal debridement: tissue sparing surgical treatment of chronic wounds. *J Wound Technol.* 2009;5:10-11.
- 7 . Jeffery SLA. Device related tangential excision in burns. *Injury, Int J Care Injured.* 2007;38S:S35-S38.
- 8. Paola LD, Brocco E, Senesi A, De Vido D, Merico M, Ninkovic S. The use of VERSAJET in the limb salvage following failure of minor amputation in diabetic foot. •Data on file report 4649/5025.
- 9. Caputo WJ, Beggs DJ, DeFede JL, Simm L, Dharma H. A prospective randomised controlled clinical trial comparing hydrosurgery debridement with conventional surgical debridement in lower extremity ulcers. *Int Wound J.* 2008;5:288-294.
- 10. Gravante G, Delogu D, *et al.* VERSAJET hydrosurgery versus classic escharectomy for burn debridement; a prospective randomized trial. *J Burn Care Res* 2007;28(5);1-5 (data for VERSAJET & escharectomy combined)
- 11. Gurunloughlu R. Experiences with waterjet hydrosurgery in wound debridement. *World J Emerg Surg* 2007 May 2;2:10

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