

Management of Advanced Skin Aging: Addressing Loss of Tone and Wrinkles with JetPeel by TavTech and JetCare Med by TavTech Serums

Performed by Laura Masini at Hautärzte Paasch Clinic, Germany

Case Description

Patient Profile: Female, Age 72.

This case demonstrates the clinical efficacy of JetPeel by TavTech together with JetCare Med by TavTech serums in treating advanced visible skin aging. The subject, a 72-year-old female, presented with deep static and dynamic wrinkles, pronounced skin laxity, uneven texture, and irregular pigmentation together representing a severe form of cutaneous aging.

During the consultation, the patient showed:

- **Marked wrinkles** (deep static and dynamic) across the forehead and periorbital areas.
- **Severe cutaneous laxity with loss of firmness and structural support.**
- **Irregular skin texture** and **uneven skin tone** consistent with advanced visible skin aging.

Treatment Process:

This protocol was designed to invigorate collagen and elastin with JetCare Med serums, targeting the core mechanisms of advanced skin aging. The treatment delivers profound improvements in skin density, elasticity and firmness.

Over the course of **two sessions**, with the second treatment performed 7 days after the first, the skin undergoes a visible transformation, regaining resilience, smoothness, and youthful vitality. The treatments ensured optimal efficacy, with no downtime and maximum patient comfort.

JetPeel by TavTech products used in JetPeel treatments

	LYMPHATIC MASSAGE	STEP 01	EXFOLIATION	STEP 02	INFUSION	STEP 03	BOOSTER INFUSION	STEP 04
TREATMENT 1	18ml JetCare Hydro		5ml JetCare Renewal Care Glycolic Acid 16%		5ml syringe of JetCare Med Relax Infusion		5ml syringe of JetCare Med Refine Boost Infusion	
TREATMENT 2	18ml JetCare Hydro		5ml JetCare Renewal Care Glycolic Acid 10%		5ml syringe of JetCare Med Refill Infusion		5ml syringe of JetCare Med Refine Boost Infusion	

TO: BEFORE JETPEEL
TREATMENTS



AFTER T2: 7 DAYS
AFTER T1



Summary and Conclusion

The treatment protocol demonstrated a powerful ability to address key skin concerns, delivering visible improvements in texture, wrinkle depth, and overall skin quality. After two sessions, the skin appeared smoother, denser, and more radiant, with noticeable gains in elasticity and firmness.

For optimal results, a series of three additional sessions is recommended, performed 7–10 days apart. This cycle maximizes the benefits by invigorating collagen and elastin, thereby improving skin density, elasticity, and firmness while enhancing overall skin resilience.

To maintain these long-term results, a maintenance plan of one treatment every 6–8 weeks is advised. This regimen sustains hydration, preserves elasticity and firmness, counteracts the appearance of wrinkles, and helps delay the deeper signs of aging — ensuring lasting skin clarity, vitality and radiance.