



THE RIGHT
WAY TO
CLEAN
AND DEBRIDE
WOUNDS

JETOX™
by TavTech




THE RIGHT WAY TO CLEAN AND DEBRIDE WOUNDS

Jetox™ by TavTech is an advanced wound cleansing, irrigation and debridement system based on world-first jet stream technology that enables thorough and painless deep cleansing treatment. Used in medical settings worldwide, the Jetox Jet lavage technology is proven to effectively cleanse complex wounds and remove non-viable tissue, optimizing the healing process and supporting medical practitioners in achieving the best possible results, with maximum comfort for patients.

JETOX™ TECHNOLOGY

CLEANING AND DEBRIDEMENT WITHOUT THE PAIN

Jetox™ by TavTech is an advanced wound cleansing, irrigation and debridement system based on world-first jet stream technology that enables thorough and painless deep cleansing treatment. Used in medical settings worldwide, the Jetox Jet lavage technology is proven to effectively cleanse complex wounds and remove non-viable tissue, optimizing the healing process and supporting medical practitioners in achieving the best possible results, with maximum comfort for patients.



CLINICAL APPLICATIONS

Jetox is designed for effective cleaning of diabetic ulcers, pressure ulcers, venous stasis ulcers, burns, and hard-to-heal tissue. The compactness, and minimal amount of solution during treatment makes it ideal for bedside application in the hospital, clinic or home settings.

IDEAL FOR USE AT:



Wound care
centers



Emergency
departments



Burn &
trauma units



Home health
services



Inpatient &
outpatient clinics



Physical therapy
departments



Podiatric
centers



Long-term
care facilities

JETOX™ SOLUTION EFFICIENCY AND EFFECTIVENESS



Thorough, efficient wound cleansing and debridement:

Promotes faster healing and enables the provision of more effective medical care.



Maximum comfort for patients:

Unique proprietary jet stream technology ensures maximum wound cleansing with minimal discomfort to the patient.



Easy, clean, compact:

Achieve outstanding results with a portable, compact device enabling easy-to-control impact. Fast simple setup and minimal cleanup required.

JETOX™ SYSTEMS ENHANCE THE LEVEL OF CARE

The Jetox™ range includes three handpieces to suit a variety of medical needs:



JETOX™ ND23

Wound cleansing system administers a gentle stream for wound cleansing and removal of non-viable tissue. The generated air creates a desensitizing effect, making it the treatment of choice for patients who are very sensitive to debridement. The system includes a flexible silicone cap that acts as a splash protector to reduce contamination.



JETOX™ ND

The mono-use disposable hydro debridement system allows a focused stream for cleansing of complex wounds and removal of biofilm while helping to decrease the bacterial load. The system includes a flexible silicone cap that acts as a splash protector to reduce contamination.



JETOX™ HDC

The Jet Lavage system for cleansing and debridement of hard to debride wounds, while simultaneously evacuating debris. The system includes built-in suction elements that easily connect to standard vacuum outlets.

TavTech Ltd is the world leader in advanced jet technology for skincare treatments and applications in aesthetic and medical care. Based on over two decades of research & development.



TAVTECH LTD
T: +972.3.536.0515
E: jetox@tav-tech.com

40 Ha'atzmaut St.
Yehud 5630417
Israel

  
www.tav-tech.com